



For more information  
on this product please  
scan QR code

---

## **Instructions For Use**

---

**Návod k použití • Gebrauchsanweisung • Instrucciones de uso  
Mode d'emploi • Οδηγίες χρήσης • Istruzioni per l'uso  
Gebruiksaanwijzing • Instruksjoner for bruk • Instrukcja użycia  
Instruções de uso • Användningsinstruktioner**

## **Disposable Oxygen Tubing**

---

**Caution - Read all instructions prior to use****Indications (Intended use)**

The Disposable Oxygen Tubing (DOT) is intended to connect the oxygen source with the oxygenation device during administration of oxygen.

**Precautions (Warnings)**

This device has been designed for single use only. Reuse / reprocessing may compromise the structural integrity, essential material or design characteristics that are critical to the overall performance of the device. This usage may lead to device failure, which may result in injury to the patient. Reuse / reprocessing or repackaging may also create a risk of contamination of the device and / or cause patient infection or cross infection. Contamination of the device may lead to injury. Do not use if package is compromised/damaged. Suitably trained healthcare personnel may only use this device. This device must be disposed of in accordance with local environmental laws and regulations.

Store 4°C to 42°C in a clean, dry environment out of direct sunlight until use.

Ensure all connections are secure and air circulates freely throughout the tube before use. Do not use oxygen mask in the proximity of any diathermy/electrosurgical or spark flame producing procedures.

Exposure to oxygen may result in explosion or fire/burns.

This device is sterile (ETO). Please note, this IFU pertains to Sterile devices sold by Fairmont Medical. If the device is non-sterile, please refer to 3rd party manufacturer for sterility specifications.

**Instructions for use**

1. Remove Oxygen Tubing from the packaging
2. Unwind the tubing
3. Attach the Oxygen tubing to the Oxygen source
4. Attach the opposite end of the tubing to the device being used (Mask)
5. Set prescribed flow on Oxygen source
6. Check Oxygen flow through the device
7. Apply device to the Patient.

**Materials**

The DOT range is made of Phthalate free PVC.  
All components and packaging are Latex free.

---



**STERILE EO**

Single sterile barrier system  
with protective packaging inside

**MD**

Is a medical device



Do not resterilise



Made in Australia



Latex free



Keep away from sunlight



Keep dry

**CE**<sup>\*</sup>  
**2460**

Notified Body



Do not use if packaging is  
damaged



Phthalate and DEHP free

**CE**<sup>\*</sup>

Conformité Européenne

**EC REP**



Fairmont Medical Products Pty Ltd  
90 Prosperity Way  
Dandenong South  
Victoria 3175  
Australia  
Ph +61 3 9720 8840  
[www.fairmontmedical.com.au](http://www.fairmontmedical.com.au)

Emergo Europe  
Westervoortsedijk 60  
6827 AT, Arnhem  
The Netherlands

PM040      Version 4.1  
Issued by    L. Wood 03/06/2024

Fairmont Medical Products Pty Ltd is certified to ISO 13485: 2016

Product specifications subject to change without notice. Pictures shown may vary from actual product.

© 2024 Fairmont Medical Products Pty Ltd all rights reserved.

\*Please note the CE 2460 mark applies to sterile versions of the product and the CE mark applies to non-sterile versions.